***Assembly and Installation Procedures for "Novasip"***

***Filter Assemblies for model numbers C\*\*\*\*P1***

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***Japanese***

***Text Required***

# **Assembly and Installation Procedures for "Novasip"**

## **Filter Assemblies for model numbers C\*\*\*\*P1**

The following procedures are required to install 'Pall' 'Novasip' disposable filter assemblies and should be thoroughly read as they contain valuable information gained by extensive experience. It is very important that all instructions are carefully followed and where appropriate be incorporated into the operator's standard operating procedures. If some of the procedures do not suit your needs, please consult Pall or your local distributor before finalising your system.

### **1. SPECIFICATIONS**

#### **Materials of construction**

Filter shell: Polyetherimide

Other materials will depend on the type of Novasip filter selected. Refer to Pall Publications USD1647 and USD1648 for further information.

#### **For use in gases and liquids**

Maximum working pressure at 40°C

in compatible fluids: 6.5 barg

Maximum differential pressure: Refer to appropriate Pall data sheet

Inlet/outlet: 1-inch/1 1/2-inch 'Triclover' compatible

Vent/drain;

Vent valve: 'Stäubli' connection to suit a RBE 03.1906 female Stäubli socket

Drain valve: Hosebarb connection to suit flexible pipe with ID between 5 and 7 mm

**WARNING: NOVASIP FILTER ASSEMBLIES SHOULD NOT BE USED WITH FLUIDS INCOMPATIBLE WITH THE MATERIALS OF CONSTRUCTION, INCLUDING CLEANING AGENTS. INCOMPATIBLE MATERIALS ARE THOSE WHICH CHEMICALLY ATTACK, SOFTEN, STRESS CRACK OR ADVERSELY AFFECT THE MATERIALS OF CONSTRUCTION IN ANY WAY. USERS SHOULD CHECK THE COMPATIBILITY OF PROCESS FLUIDS BEFORE USE.**

It is the user's responsibility to check actual operating conditions to ensure that the filter assembly is compatible with the application and within local safety codes.

## **2. RECEIPT OF EQUIPMENT**

- (a) Store the filter assemblies in clean, dry conditions, out of direct sunlight, and wherever practical in the packaging as delivered.
- (b) DO NOT remove from its protective packaging until just before installation.
- (c) Ensure that the type of the filter assembly selected is suitable for the application. In addition to the part number each filter assembly is identified by an unique identification batch/serial number.

**IMPORTANT: NOVASIP FILTER ASSEMBLIES SHOULD BE INSPECTED FOR SIGNS OF DAMAGE BEFORE USE. IF USED OVER EXTENDED TIME PERIODS, IT IS RECOMMENDED THAT INSPECTION IS PERFORMED AT LEAST EVERY WEEK. IF DAMAGE OR ANY OTHER POTENTIAL MALFUNCTION IS OBSERVED, THE FILTER ASSEMBLY SHOULD BE REPLACED.**

## **3. INSTALLATION AND OPERATION**

### **3.1 Installation**

Before installation it is essential to verify that the filter type selected is suitable for the gas or liquid to be filtered and to follow the appropriate instructions listed below.

**Important: In liquid and pressurised gas applications the filter assembly must be installed with the process flow as indicated by the flow direction arrow on the assembly.**

- (a) Install the filter assembly in line using appropriate connections.
- (b) The filter assembly should be positioned in an attitude that allows effective venting and draining of the filter to be carried out before and during operation.
- (c) The filter assembly should be installed in an appropriate orientation allowing integrity testing as required.
- (d) Where a positive pressure exists downstream of the filter, a sensitive check valve may be needed to prevent back pressure damage due to reverse flow.
- (e) Where a pulsating flow is present, the filter should be protected by a surge tank upstream.
- (f) Where a rapidly closing downstream valve is present, the possibility of pressure pulsing and subsequent filter damage exists. The filter should be protected by a surge tank between valve and filter.
- (g) Side and end loads should be avoided during installation and use.
- (h) Allowance should be made for expansion during sterilisation.

- (i) Overtightening of the inlet and outlet Triclover clamps may result in damage of the inlet and outlet Triclover connectors at steaming temperatures. It is recommended that the clamps are fully tightened by hand, and then loosened by one turn. It is also recommended that users verify that this provides a leak proof seal.

**WARNING: NOVASIP FILTER ASSEMBLIES HAVE BEEN EXTENSIVELY VALIDATED FOR USE IN PRESSURISED GAS SYSTEMS AND FOR STEAM STERILISATION IN PLACE. USERS SHOULD TAKE THE APPROPRIATE PRECAUTIONS ASSOCIATED WITH SUCH PRESSURISED AND HIGH TEMPERATURE SYSTEMS TO PROTECT OPERATORS SUCH AS SAFETY GLASSES AND GLOVES. IN ADDITION PALL RECOMMEND THE USE OF A PROTECTIVE SHIELD TO PROTECT OPERATORS FROM LIVE STEAM IN THE EVENT OF A LEAK OR BREAKAGE.**

### **3.2 Operation**

On installation and prior to steaming verify the integrity of the assembly.

In liquid applications the filter assembly must be vented by opening the vent valve when the filter is initially brought on stream. The vent valve should then be closed when liquid issues from it.

## **4. INTEGRITY TESTING**

Filter integrity should be verified before and after filtration by means of an appropriate integrity test procedure. Depending on the filter type and grade, the Pall Forward Flow, water intrusion and bubble point tests may be suitable.

### **4.1 Forward Flow Test**

The Pall Forward Flow test can be performed on a sterile or non-sterile filter assembly. The test is suitable for all types of Pall sterilising grade and other fine particulate grade filters.

In summary the test requires establishment of a pre-determined test pressure using compressed air or nitrogen in the filter assembly, on the upstream side of the pre-wetted filter membrane. Under the influence of this constant pressure, there is a flow of test gas through the membrane resulting from diffusive and open pore flow. Pall has correlated measurement of the volume of gas flow with the results of bacterial retention tests. The validated integrity test enables filter performance capability to be determined in a pharmaceutical process without the risk of contaminating the test system. Pall filter validation guides carry more detail for specific filter types.

The 'Palltronic' 'Flowstar' or 'TruFlow' test instruments carry out the Forward Flow test by measurement of diffusive flow at constant pressure from the upstream side of the filter, thereby eliminating the risk of contaminating the downstream (sterile) side of the filter. The Palltronic instrument is connected to the filter assembly after wetting the filter with an appropriate test liquid. After entering the necessary information and test parameters the instrument will carry out the Forward Flow test and report the result. During the test temperature should not vary by more than  $\pm 1^{\circ}\text{C}$ .

Appropriate wetting liquids for the Forward Flow test with Pall filters are as follows:

Filter Medium Type	Wetting Liquids
Ultipor N <sub>66</sub>	Water, 60/40 (V/V) IPA*/water
N <sub>66</sub> Posidyne	Water, 60/40 (V/V) IPA*/water
Fluorodyne II	Water, 60/40 (V/V) IPA*/water
Ultipor VF DV50	30/70 (V/V) IPA*/water
Emflon PFR	60/40 (V/V) IPA*/water, 25/75 t-butyl alcohol water

\* Isopropyl alcohol (IPA)

During process operation it is frequently more convenient to be able to test a filter whilst wet with the actual process fluid. This is nearly always possible; Pall Scientific and Laboratory Services should be contacted for further information.

#### 4.2 Water Intrusion Test

The Pall Water Intrusion test can be performed on a sterile or non-sterile filter assembly. The test is suitable for Pall 'Emflon' PFR hydrophobic sterilising grade gas filters.

If the upstream side of a dry hydrophobic filter assembly is filled with water and pressurised, the hydrophobic nature of the porous membrane will prevent the bulk flow of liquid through the filter until the intrusion pressure is reached. At pressures below the intrusion pressure a small but measurable flow of water through the membrane occurs. The presence of larger pores in the filter will be detected by an increased flow resulting from to bulk water flow through these pores. This principle forms the basis of the water intrusion test. Pall has correlated measurement of water flow with the results of bacterial retention tests, permitting the test to be used to determine filter performance capability. Pall Scientific and Technical Reports carry more detailed information about this test.

The upstream side of the filter assembly is filled with distilled or deionised water at 20°C ± 2°C and of surface tension not less than 0.071 N/m. Under the influence of a pre-determined applied pressure, there is a fall in water level in the filter assembly due to factors such as pleat compression and passage of trapped air through the membrane. After allowing an appropriate period for stabilisation, water flow is measured while maintaining constant upstream pressure.

The Palltronic Flowstar or TruFlow test instruments carry out the water intrusion test by measurement of flow from the upstream side of the filter, thereby eliminating the risk of contaminating the downstream (sterile) side of the filter. The Palltronic instrument is connected to the filter assembly after filling the filter assembly with water. After entering the necessary information and test parameters the instrument will carry out the water intrusion test and report the result. On completion of the test the filter assembly is drained of water and, in many applications, the filter is now ready for use. A drying procedure can be used in those applications where the filter must be dry before use.

### 4.3 Bubble Point Test

A bubble point type test can be performed on a sterile or non-sterile filter assembly. The test can be applied to all types of Pall sterilising grade and other fine particulate grade filters.

A progressively increasing pressure is applied using compressed air or nitrogen in the filter assembly, on the upstream side of the pre-wetted filter membrane. A flow of gas occurs through the membrane resulting from expulsion of the wetting liquid from the largest filter pores.

The Palltronic Flowstar or TruFlow test instruments carry out a bubble point test from the upstream side of the filter, thereby eliminating the risk of contaminating the downstream (sterile) side of the filter. The Palltronic instrument is connected to the filter assembly after wetting the filter with an appropriate test liquid. After entering the necessary information and test parameters the instrument will carry out a bubble point test and report the result.

Appropriate wetting liquids for a bubble point test with Pall filters are as follows:

Filter Medium Type	Wetting Liquids
Ultipor N <sub>66</sub>	Water, 60/40 (V/V) IPA*/water
N <sub>66</sub> Posidyne	Water, 60/40 (V/V) IPA*/water
Fluorodyne II	Water, 60/40 (V/V) IPA*/water
Ultipor VF DV50	30/70 (V/V) IPA*/water
Emflon PFR	60/40 (V/V) IPA*/water, 25/75 t-butyl alcohol water

During process operation it is frequently more convenient to be able to test a filter whilst wet with the actual process fluid. This is nearly always possible; Pall Scientific and Laboratory Services should be contacted for further information.

### 5. STERILISATION

THESE DISPOSABLE NOVASIP FILTER ASSEMBLIES MAY BE STEAM STERILISED IN PLACE OR BY AUTOCLAVE

**WARNING: DO NOT AUTOCLAVE FILTER ASSEMBLY IN THE BAG.**

5.1 Novasip filter assemblies are supplied non-sterilised.

5.2 The Novasip filter assembly should only be steam sterilised or autoclaved as follows:


- (a) Up to 142°C
- (b) DO NOT steam sterilise at higher temperatures.
- (c) For specific steam sterilisation limits please contact Pall.
- (d) Quench cooling is not recommended

5.3 Detailed information on sterilisation of Pall filters is provided in Pall Publication USD805j and a typical procedure for steam in place sterilisation of a dry filter is given below. Further information can be obtained by contacting Pall.

5.4 For gas filter assemblies a numbered plastic ring is supplied, for fitting if required, this can be used to record the number of steam sterilisation cycles carried out.

*Typical Procedure for Steam in Place Sterilisation of a Dry Filter Assembly*

A recommended filter installation is presented in Figure 1. Pressure gauges which can be read with accuracy over the range of 0-3 bar must be installed to monitor steam pressure and differential pressure across the filter assembly during the sterilisation cycle. To ensure effective sterilisation, steam pressure (measured at position T) in the assembly should be held at a minimum of 121°C (1.1 bar g saturated steam) for the minimum time validated by the user as necessary to achieve system sterilisation.

*Figure 1. Recommended filter installation for steam in place sterilisation of a dry filter assembly*

1. Follow installations instructions in section 3.1 and ensure all valves are closed.
2. Fully open valve C.
3. Fully open condensate drain trap or valve I, housing drain valve J and housing vent valve G.
4. Preset steam pressure (P4) to 300 mbar above the steam pressure required at the filter assembly. After condensate has been expelled from I partially close valve I (if necessary).
5. Slowly open steam valve B to admit steam to system. After condensate has been expelled from J, partially close valve J.
6. Partially close vent valve G when steam flow is evident, ENSURING THAT PRESSURE AT P2 REMAINS WITHIN 300 MBAR OF PRESSURE AT P1. Partially open drain valve D to drain condensate.
7. Permit steam to flow through the system until steam pressure is stabilised, and adjust the regulated steam supply until the validated temperature is achieved at position T. Monitor temperature at T for the necessary sterilisation time.

**ENSURE THAT PRESSURE AT P2 REMAINS WITHIN 300 MBAR OF PRESSURE AT P1. IT IS RECOMMENDED THAT STEAM STERILISATION IS FOLLOWED BY AIR BALLASTING AS DETAILED IN THE FOLLOWING SECTION.**

8. Preset pressure (P3) of regulated air or N2 at 200 mbar above steam pressure (P4).
9. When sterilisation time is complete close drain valves D, J and I and vent valve G. Close steam valve B and immediately introduce preregulated air or N2 through valve F.
10. To assist cooling steam may be flushed from the assembly by carefully opening vent valve G and drain valve J. Close valves G and J after flushing.
11. Allow assembly to cool to ambient or to process fluid temperature.
12. Close air or nitrogen valve F.
13. Relieve the gas pressure in the filter assembly via vent valve G. Filter assembly is now ready for use.

## **6. FILTER ASSEMBLY REPLACEMENT**

A filter assembly should be replaced when the maximum allowable differential pressure of the filter has been reached, (refer to appropriate Pall data sheet), the flow rate has become unacceptable, the steam life of the filter has been exceeded or damage is evident, whichever occurs first. Discard filter assembly in accordance with local Health and Safety Procedures. No attempt should be made to clean disposable filter assemblies.

## **7. SCIENTIFIC AND LABORATORY SERVICES**

Pall operate a technical service to assist in the application of all filter products. This service is available to you without charge and we welcome your questions so that we can help. In addition, a full network of technical representatives are available throughout the world.

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# Japanese Films

# Customer support throughout the world

## Pall has sales offices in the following countries:

**Australia:** Pall Australia, Adelaide, Brisbane, Mackay, Melbourne, Newcastle, Perth and Sydney.

**Austria:** Pall Austria, Vienna. **Belgium:** Pall bvba., Machelen (Brussels).

**Brazil:** Pall Industrial do Brasil Ltda., Sao Paulo. **Canada:** Pall Canada Ltd., Ontario.

**China:** Pall Filter (Beijing) Co. Ltd., Beijing.

**France:** Pall Industrie, Paris. **Germany:** Pall GmbH Filtrationstechnik, Frankfurt.

**Hong Kong:** Pall Fluid Pte Ltd., Wanchai. **India:** Pall Pharmalab Pvt. Limited, Bombay. **Italy:** Pall Italia s.r.l., Milan.

**Japan:** Nihon Pall Ltd., Tokyo. **Korea:** Pall Korea Ltd., Seoul. **Poland:** Pall Poland Ltd., Warsaw.

**Russia:** Pall Filtrationstechnik GmbH, Moscow. **Singapore:** Pall Filtration Pte Ltd. **Spain:** Pall España s.a., Madrid.

**Switzerland:** Pall (Schweiz) AG, Muttenz. **Taiwan:** Pall Filtration Pte Ltd., Taipei.

**United Kingdom:** Pall Ultrafine Group, Portsmouth and Glasgow. **USA:** Pall Corporation, New York.

## Pall BioPharmaceuticals is also represented by distributors in:

**Croatia:** Ipros, Ljubljana. **Denmark:** Colly Company a/s, Copenhagen. **Egypt:** Akhnaton Trading & Representation, Cairo.

**Finland:** Oy Colly Company ab, Vantaa. **Greece:** Th. Chalimas EPE, Athens. **Iran:** Panaid Engineering Limited, Tehran.

**Ireland:** Gilroy Filtration, Belfast/Cork/Dublin. **Israel:** B.T.S. Baruch & Company Technologies & Supplies Limited, Herzliya.

**Macedonia:** Alfaprodukt, Skopje. **Netherlands:** Hitma B.V., Uithoorn.

**New Zealand:** A.M.E. Industries (NZ) Limited, Hamilton. **Norway:** Caspar Hansen A/S, Hovik.

**Pakistan:** Kiran Enterprises (Pvt) Limited, Karachi. **Portugal:** Semortrade - Importacao E Exportacao Lda, Lisbon.

**Saudi Arabia:** Saudi Arabian Industrial & Trading Establishment, Dammam. **Slovenia:** Ipros, Ljubljana.

**South Africa:** Deomed (Pty) Limited, Johannesburg. **Sweden:** Colly Company ab, Stockholm.

These Accreditations/Approvals have been awarded to Pall Europe Limited

Pall Europe Corporate Services Ltd.  
Scientific & Laboratory Services

QUALITY MANAGEMENT SYSTEMS

871326 921229 871806 921227 921228 923031  
871805 921231 871807 921233

ENVIRONMENTAL MANAGEMENT SYSTEMS

770048

## BioPharmaceuticals

### Pall BioPharmaceuticals

A Division of Pall Europe Limited  
Europa House,  
Havant Street,  
Portsmouth PO1 3PD England

+44(0)1705 303303 phone  
+44(0)1705 302506 fax

Visit us on the web at [www.pall.com](http://www.pall.com)

e-mail: [UltrafineUK@pall.com](mailto:UltrafineUK@pall.com)

Because of developments in technology these data or procedures may be subject to change. Consequently we advise users to review their continuing validity annually.