

Biotech

Validation Guide

USTR 3393

# iCELLis<sup>®</sup> 500+ GenQ Bioreactors

# 1 Introduction

## 1.1 iCELLis 500+ GenQ Bioreactor

This validation guide contains data applicable to iCELLis 500+ GenQ bioreactors. An iCELLis 500+ GenQ bioreactor is an automated single-use (SU) bioreactor that provides excellent cell growth conditions for adherent cells. Central to the iCELLis 500+ GenQ bioreactor technology is the use of a compact fixed bed, filled with custom macrocarriers.

Evenly distributed media circulation is achieved by a built-in magnetic drive impeller, ensuring low shear stress and high cell viability. The cell culture medium flows through the fixed bed from the bottom to the top. At the top, the medium falls as a thin film down the outer wall where it takes up  $O_2$  to maintain high  $k_{L}a$  in the bioreactor. This unique waterfall oxygenation, combined with a gentle agitation and biomass immobilization, enables the compact iCELLis 500+ GenQ bioreactor to achieve and maintain high cell densities, equalling the productivity of much larger stirred tank units.

The iCELLis 500+ bioreactor system consists of two major components:

- a control system
- a single-use bioreactor

This guide focuses on the second element, which is a gamma irradiated plastic bioreactor that provides cells with a growth area up to 500 m<sup>2</sup>. This document describes testing demonstrating the suitability of iCELLis 500+ GenQ bioreactors to culture cells in the biotechnology industry.

The single-use iCELLis 500+ GenQ bioreactor is composed of the following elements:

- 2 probe supports for temperature (T), dissolved oxygen (DO) and pH sensors
- Lid containing 8 connections:
  - o 5 inlets: gas inlet, 2 feed inlets, inoculum inlet, base inlet
  - o 3 outlets: gas outlet, feed outlet, harvest/drain outlet
- A sampling line (allows to sample before and after fixed bed)
- A fixed bed support integrating a pump housing, filled with macrocarriers and a top grid to retain the macrocarriers
- An integrated magnetic-driven impeller and a pump housing ring
- An outer shell including a double jacket at the bottom for temperature regulation
- A biomass probe (optional)

The purpose of this technical note is to provide specific instructions that should be followed while using the iCELLis 500+ GENQ single-use bioreactors with Hamilton\* probes with the iCELLis 500+ control system.



In this document, "iCELLis 500+ GenQ bioreactor" refers to iCELLis 500+ Q generation single-use bioreactors. Please refers to Customer Notifications CN-RPII-BJGPN9 and CN-RPII-BLWMAK for more information.

PALL Biotech

2

### Figure 1

The single-use iCELLis 500+ GenQ bioreactor



The bioreactor has two main product contact materials:

- Polyethylene terephthalate (PET) used for the macrocarriers (98% of the total product contact surface area for the 500 m<sup>2</sup> version, 85% for the 66 m<sup>2</sup> version)
- Polyethylene terephthalate glycol (PETG) used for the main housing components (1% of the total product contact surface for the 500 m<sup>2</sup> version; 6% for the 66 m<sup>2</sup> version)
- The same manufacturing methodology and same component materials are used to produce different types of iCELLis 500+ GenQ bioreactors (66, 100, 133, 200, 333 or 500 m<sup>2</sup> carriers, with or without biomass probe). Therefore, the validation results mentioned in this guide are applicable for all container models.

### 1.2 History

In order to meet customer requests and to improve the performance of the bioreactor, changes were made to the existing disposable bioreactor and a new Q generation of iCELLis 500+ single-use bioreactors was launched.

The changes include:

- Complete transition from Polestar to Hamilton\* pH and DO sensors
- Addition of a sampling line
- Change to injection molded fitments
- Change of the macrocarrier shape

# 2 Material Compliance

## 2.1 BSE/TSE

All product contact materials of the iCELLis 500+ GenQ bioreactor components are certified ADCF (Animal Derived Component Free) or comply with the requirements of EMA 410/01 rev 3- July 2011.

# 2.2 Biological Safety Tests

The purpose of these tests is to evaluate the biological suitability of the construction materials of iCELLis 500+ GenQ bioreactors. All polymer-based product contact materials of the iCELLis 500+ GenQ bioreactor are certified USP <88> Class VI. The main product contact materials of the bioreactor (PET and PETG) are also certified USP <87> and USP <661>.

# 3 Extractables

Extractables are chemical entities that migrate from process equipment under appropriate exaggerated conditions (e.g. solvent, temperature) that exceed worst case process conditions. They represent chemicals that may (or may not) be found in actual process fluid, final formulations and/or finished dosage.

Polymers, used in the production of pharmaceutical containers and medical devices, cannot be considered as pure compounds. Medical grade polymers should be a blend of base polymers with a broad range of chemicals that may be added for several reasons.

Polymer additives are added to improve the processability of the polymer to enhance its end-use performance in various ways. Moreover, residues of unexpected and undesirable compounds may also be present and can affect the biocompatibility and the toxicity of the materials. Residues can originate from unreacted monomers, solvents used in the process, polymerization catalysts, surfactants or polymer degradation products.

Both the pharmaceutical and medical device industries require to assess the toxicological risk associated with the use of the product. The FDA as well as the EMEA suggest in these cases to perform an extraction study on pharmaceutical containers or medical devices to determine which chemical species may migrate out of the polymer material and at what concentration level. These data will allow evaluation of the potential toxicological risk cause by the extracted substances.

# 3.1 Selection of Test Materials

When defining the extractables program for the iCELLis 500+ GenQ bioreactor we started with a materials risk evaluation of the product contact materials. From a risk evaluation point of view, the most important materials are those with the largest contact surface, largest contact time and those which might represent a safety risk. This approach reduces the list of materials to be tested to the following 2 materials:

- Polyethylene terephthalate glycol (PETG) main contact surface in the bioreactor (sum of vessel, pump ring, pump housing basket support, fixed bed plates and lid).
- Polyethylene terephthalate (PET) main contact surface in the bioreactor.

For the other product contact materials, the rationale for not testing their extractables is due to their small contact surface and/or short contact time.

# 3.2 **Pre-Treatment Prior to Extraction**

Gamma irradiation can potentially increase the amount of extractables, so the test articles should be exposed to this sterilization treatment at the maximum dosage allowed prior to extraction. For this reason, test articles were gamma irradiated at 50 kGy prior to extraction.

## 3.3 Extraction Conditions

Extractions should be performed at exaggerated conditions, considered to be harsher than normal process conditions. To prevent over-prediction of extractable compounds however, conditions should not be overly extreme.

#### Table 1

Extraction conditions

Solvent	Temperature and Time
10% ethanol in water*	8 hour, reflux

\* This solvent was chosen because because the iCELLis 500+ GenQ bioreactor is used in applications involving aqueous solutions. Using stronger solvents will generate a far too exaggerated condition.

### 3.4 Analysis of the Extract

To identify the extract, the following analytical methods were applied:

Headspace Gas Chromatography/Mass Spectrometry (HS-GC/MS) to determine volatile organic components (VOC)

Volatile organic molecules which migrate into the contact solution during a prolonged contact step between the test material and the extraction solution, may come from various sources, for example:

- monomer residues
- solvent residues from various production or polymer treatment steps (e.g. washing)
- smaller polymer breakdown products

The selected technique for this analytical method – Purge and Trap coupled to a Gas Chromatography (GC) with Mass Spectrometry (MS) as a detection technique – allows to carry out an identification of the target analytes based on both the retention time of the analytes in the chromatogram and the mass spectrum of the eluting compound at this specific retention time. The concentration of detected volatile compounds could be estimated by a semi-quantitative internal calibration method. For this purpose, the relative analytical response of a reported compound was compared to the response of an appropriate internal standard toluene-d<sub>8</sub>, added to the samples at a fixed concentration.

# Gas Chromatography/Mass Spectrometry (GC/MS) to determine semi-volatile organic components (SVOC)

A lot of potential organic migration products are not volatile enough to be detected via Purge and Trap GC/MS. These thermostable compounds can be studied via GC/MS, even though under different experimental and instrument condition than used with Purge and Trap GC/MS. These products are called the semi volatile compounds and may come from various sources, for example:

- process lubricants
- plasticizers
- antioxidants
- polymer degradation products
- solvents with a higher boiling point

Dichloromethane extraction followed by GC with MS can identify the target analytes based on both the retention time of the analytes in the chromatogram and the mass spectrum of the eluting compound at this specific retention time. The concentration of detected semi-volatile compounds

could be estimated by a semi-quantitative internal calibration method. The relative analytical response of a reported compound was compared to the response of an appropriate internal standard (2-Fluorobiphenyl), added to the samples at fixed concentration.

In a first step the water samples are extracted with an organic solvent with a low boiling point. In general, the solubility of organic compounds is much larger in organic solvents than in water. During the extraction, most organic compounds are concentrated in the organic phase.

# *Liquid Chromatography/Mass Spectrometry (LC/MS) to determine non-volatile organic components (NVOC)*

For migration compounds that are non-volatile or non-thermostable (like antioxidants, fillers, plasticizers, polymerization or hydrogenation catalysts, anti-slip agents and other polymer additives), Liquid Chromatography (LC) is better suited as an analytical tool compared to GC. Mass spectral detection was used because it offers numerous advantages over the traditional LC/UV-VIS technique, such as additional molecular and structural information of the compounds, a higher sensitivity of the instrument, and a better identification of the target compounds. The analytes present in the extract are identified from their retention time and the corresponding mass spectrum in a semi-quantitative manner. The concentration of detected non-volatile compounds could be estimated by a semi-quantitative internal calibration method.

The relative analytical response of a reported compound was compared to the response of an appropriate internal standard (Tinuvin<sup>•</sup> 327). A compound specific relative response factor (RRF) is available for several target compounds. This RRF is used for the more accurate quantification of identified target compounds. For other compounds the RRF is assumed to be 1.

In a first step the water samples are extracted with an organic solvent with a low boiling point. In general, the solubility of organic compounds is much larger in organic solvents than in water. During the extraction, most organic compounds are concentrated in the organic phase.

### 3.5 Results

The migration behavior was verified and compared with the blank solution, stored in a glass bottle. In this way, it becomes possible to discriminate migration and material degradation compounds from compounds originating from the matrix/solution degradation under test conditions.



### Table 2

Extractables results for polyethylene terephthalate glycol (PETG), minimum 50 kGy gamma irradiated\*.

Analysis	Compound	Concentration (µg/g)
Volatile organic compounds (neat). Reporting limit: 0.50 µg/g.	No compounds detected above the reporting limit and different from the blank	N/A
	Unknown	0.11
	Ester of terephthalic acid	0.067
	Bis (2-hydroxyethyl) – terephthalate	0.057
Semi volatile organic compounds	PET dimer compound	0.38
Reporting limit: 0.050 µg/g.	PET dimer related compound	0.39
	Probably PET related	0.27
	Unknown	0.20
	PET dimer	0.59
	PET dimer	0.75
	Probably PET related	0.053
	Probably PET related	0.081
	PET dimer	0.089
	PET trimer	0.63
	PET related compound	0.055
Non-volatile organic compounds. Reporting limit: 0.050 µg/g. Positive ionization mode.	PET linear vinyl (CDHM) dimer of PET cyclic (CDHM glycol diglycol) dimer	0.24
	PET (CDHM) cyclic trimer	0.052
	PET cyclic trimer	0.82
	PET linear bisglycol dimer	0.81
	PET linear (monoglycol) dimer methyl ester	0.096
	Unknown	0.12
	Unknown	0.09
	PET dimer	0.12
	PET trimer	2.2
	PET linear (CDHM monoglycol) dimer	0.41
	Unknown	0.055
	Unknown	0.11
	PET trimer	0.15
	PET related compound	0.27
	PET related compound	0.056
	PET related compound	0.062
	PET (CDHM) cyclic trimer	0.39
Reporting limit: 0.050 ug/a.	Unknown	0.052
Positive ionization mode.	PET cyclic trimer	1.5

\* Surface to volume ratio: 1 g/10 mL. Results are reported in  $\mu$ g/g units, meaning the amount of extractable compound ( $\mu$ g) per amount of extracted polymer material (g).

### Table 3

Extractables results for polyethylene terephthalate (PET), non-woven, minimum 50 kGy gamma irradiated\*.

Analysis	Compound	Concentration (µg/g)
	Acetaldehyde	0.93
	Methyl formate	2.3
Volatile organic compounds (neat)	Formic acid	1.0
Reporting limit: $0.50 \ \mu g/g$ .	Acetic acid	2.2
	Diethyl terephthalate	0.23
	Ester of terephthalic acid	0.36
	Ester of terephthalic acid	0.52
	Bis(2-hydroxyethyl)-terephthalate	1.1
	PET dimer	3.4
Semi-volatile organic compounds	PET dimer	5.6
Reporting limit: 0.10 µg/g.	Cyclic PBT dimer	0.75
	Probably PET related	0.31
	Unknown	0.13
	PET dimer	0.65
	PET dimer	1.2
	PET dimer	7.4
	Probably PET related	0.86
	PET dimer	0.96
	Probably PET related	0.66
	PET dimer	0.82
	PET trimer	0.20
	PET dimer	0.10
	PET trimer	0.11
Non-volatile organic compounds.	PET cyclic trimer	4.5
	PET cyclic tetramer	0.27
Positive ionization mode.	PET cyclic pentamer	0.26
	2-hydroxyethyl terephthalate	0.13
	Impurity related to PBT	0.45
	PET dimer	2.6
	PET dimer	0.19
	PET dimer	21
	Unknown	1.4
	PET dimer	0.90
	Probably PET related	0.29
	Unknown	0.92
	Unknown	0.79
New collectile environments of the second second	PET trimer	1.2
Reporting limit: 0.10 ug/g	PBT dimer	0.16
Negative ionization mode.	PET cyclic trimer	6.9

\*Surface to volume ratio: 1 g/mL. Results are reported in  $\mu$ g/g units, meaning the amount of extractable compound ( $\mu$ g) per amount of extracted polymer material (g).



8

# 4 Cleanliness

To meet the cleanliness specifications for the operational iCELLis 500+ GenQ bioreactor, described in US and European Pharmacopoeia, Pall Biotech uses internal specifications that must be met by all components of the bioreactor.

### 4.1 Particles

In accordance with Ph. Eur. 2.9.19 and USP <788>, a light obscuration particle count was performed on rinse volume samples from an empty iCELLis 500+ GenQ bioreactor and macrocarriers to determine the particulate matter.

#### Table 4

Particle concentration in filled bioreactor\* (particles/mL)

	Particle Size	
	≥ 10 µm	≥ 25 µm
Average of 3 batches of macrocarriers	488	179
Average of 3 empty iCELLis 500+ GenQ bioreactors	17	3
Sum	505	182

\* Considering 3600 g macrocarriers in a volume of 65 L in the iCELLis bioreactor model with 500 m<sup>2</sup>

These results were generated for information purposes by Pall's R&D facility in Hoegaarden, Belgium. During production no cleanliness data is measured as a release criterium.

# 5 Sterility

The iCELLis 500+ GenQ bioreactor is not shipped with a sterility claim, but each bioreactor is shipped with an irradiation certificate for minimum 30 – maximum 50 kGy.

# 5.1 Bioburden Prior to Gamma Irradiation

The bioburden load prior to gamma irradiation of a representative iCELLis 500+ GenQ bioreactor and of the macrocarriers were analyzed in accordance with ISO11737-1.

#### Table 5

Bioburden content in bioreactor\* (CFU)

	Bioburden Content in Bioreactor* (CFU)
Average of 3 batches of macrocarriers	5810
Average of 3 empty iCELLis 500+ GenQ bioreactors	< 2817
Sum	< 8627

\* Considering 3600 g macrocarriers in a volume of 65 L in the iCELLis 500+ GenQ bioreactor model with 500 m<sup>2</sup>

These results were generated for information purposes by Pall's R&D facility in Hoegaarden, Belgium. During production no bioburden data is measured as a release criterium.

### 5.2 Dose Mapping

All iCELLis 500+ GenQ bioreactors are shipped with an irradiation certificate for minimum 30 kGy – maximum 50 kGy. According to ISO11737-2, this results in a sterility assurance level (SAL) of  $10^{-6}$  for items with an average bioburden load < 23000 CFU.

# 6 Shelf Life

An 18-month shelf life is ensured for gamma-irradiated iCELLis 500+ GenQ bioreactors based on 3 aging studies:

- A real time aging study of the plasma treatment confirms hydrophilicity of the macrocarriers to be within internal specifications after 18 months.
- A real time aging study of DO and pH sensors, performed by the supplier, which results in a shelf life of 18 months after gamma irradiation. A study is ongoing by Pall's R&D facility in Hoegaarden, Belgium to confirm this statement.
- An accelerated aging study of iCELLis 500+ GenQ bioreactor is performed in accordance with ASTM F1980-16 and demonstrates conformity with the specifications after 18 months for the following parameters:
  - Bioreactor integrity
  - Fitment integrity
  - o Sampling line quality
  - o Sensor functionality and integrity
  - PHBS connection quality

# 7 Transportation Study – Packaging

The bioreactor packaging is tested according to ISTA 3E and confirmed to be compliant. These tests include shock, compression, and vibration tests.

# 8 Sensor Accuracy

In case of a one-point calibration in the range from 6.5 to 8.5, an accuracy of  $\pm$  0.1 is guaranteed for the pH sensor after gamma irradiation and with product calibration.

For the DO sensor, an accuracy of  $\pm 3.5\%$  between 15 and 75% and of  $\pm 7\%$  between 75% and 150% after one-point calibration is provided after gamma irradiation and with product calibration.

# 9 Fixed Bed Compaction

Pall provides all bioreactors with a uniform carrier distribution within the fixed bed for optimal cell growth and cell distribution.

10 PALL Biotech



# Biotech

#### **Corporate Headquarters**

Port Washington, NY, USA +1.800.717.7255 toll free (USA) +1.516.484.5400 phone

#### European Headquarters

Fribourg, Switzerland +41 (0)26 350 53 00 phone

Asia-Pacific Headquarters Singapore +65 6389 6500 phone

#### Visit us on the Web at www.pall.com/biotech

#### Contact us at www.pall.com/contact

#### International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States, and Venezuela. Distributors in all major industrial areas of the world. To locate the Pall office or distributor nearest you, visit www. Pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

© 2020 Pail Corporation. The Pail logo, Pail and iCELLis are trademarks of Pail Corporation. © indicates a trademark registered in the USA and TM indicates a common law trademark. Filtration. Separation. Solution is a service mark of Pail Corporation. • Hamilton is a trademark of Hamilton Company and Tinuvin is a trademark of BASF SE Company