




# GeneDisc® Method for Spoilage Yeast in Beverages: A Guide to Testing Strategies

## Benefits

<b>Accelerated decision-making</b>	Enables early preventive controls to reduce cost of product scrap, product recall or additional product processing related to product spoilage. Speeds up batch release to reduce storage cost.
<b>Fast corrective actions implementation</b>	Reduces negative financial impact of spoilage once detected with rapid root cause analysis.
<b>Adaptable informative method designed for beverage industries</b>	Gets relevant information with two GeneDisc Plates (Yeast Screening and Yeast ID) and three testing strategies (analysis with enrichment, direct monitoring, and <i>Brettanomyces</i> quantification).
<b>Reduced hands-on cost</b>	Ease of use and yeast identification information simplifies testing workflows and on-site implementation.

## Assess Contamination in 2 hours




### Technical Information

Sensitivity	Filterable samples: As low as 1 cell/mL Unfilterable samples: As low as 85 cells/mL
Time to Results	Reduced to 2 hours
Plate Options	<ul style="list-style-type: none"> <li>• Yeast Screening</li> <li>• Yeast ID for identification of the 12 major spoilage yeast genera and species simultaneously</li> </ul>
Internal Positive Control	To ensure result accuracy, each sample analysis includes an internal positive control.

**When quick results are your priority**

## Reach High Sensitivity



### Technical Information

Sensitivity	Down to 1 cell / sample
Enrichment	As low as 28 hours
Time to Results	Enrichment time + 2 hours
Plate Options	<ul style="list-style-type: none"> <li>• Yeast Screening</li> <li>• Yeast ID for identification of the 12 major spoilage yeast genera and species simultaneously</li> </ul>
Internal Positive Control	To ensure PCR result accuracy, each sample analysis includes an internal positive control.

**When precise information is your priority**

## Monitor *Brettanomyces* Level in 2 hours


Cell  
concentration



Cell  
lysis



PCR  
analysis



### Technical Information

Quantified Targets	<i>Brettanomyces</i> spp. and <i>Brettanomyces bruxellensis</i>
Sample Types	Designed for wine process samples from grape must to bottling
Sensitivity	As low as 1 cell /mL
Quantification Range	As low as 1 to 100,000 cells/mL
Time to Results	Reduced to 2 hours
Detected Targets	Allows simultaneous detection of 10 additional spoilage yeast genera and species
Internal Positive Control	To ensure result accuracy, each sample analysis includes an internal positive control.

**To preserve wine sensory characteristics**



Pall Corporation

Pall GeneDisc Technologies

25 Harbor Park Drive  
Port Washington, NY 11050  
+1 516 484 3600 telephone  
+1 866 905 7255 toll free US


genedisc@pall.com

Visit us on the Web at [www.pall.com/genedisc](http://www.pall.com/genedisc)

#### International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States, and Venezuela. Distributors in all major industrial areas of the world. To locate the Pall office or distributor nearest you, visit [www.pall.com/contact](http://www.pall.com/contact).

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

© 2017, Pall Corporation. Pall, , and GeneDisc are trademarks of Pall Corporation. ® indicates a trademark registered in the USA and TM indicates a common law trademark. *Filtration.Separation.Solution.sm* is a service mark of Pall Corporation.