## Acrodisc® Syringe Filters

For HPLC Sample Preparation and Solvent or Aqueous Filtration

### HPLC Certification

Pall Life Sciences certifies the following HPLC Acrodisc syringe filters have been tested for compatibility with common HPLC solvents using established HPLC procedures. To verify low levels of UV-detectable extractables, representative samples are tested using a highly sensitive HPLC technique.

**Note:** Extractable levels may be reduced by flushing the filter prior to use.

### Ordering Information

Part Number	Description	Pkg PK/Case
4554	GHP Acrodisc Minispike, 0.2 µm, 13 mm	100/300
4556	GHP Acrodisc Minispike, 0.45 µm, 13 mm	100/300
4564	GHP Acrodisc, 0.2 µm, 25 mm	50/200
4560	GHP Acrodisc, 0.45 µm, 25 mm	50/200
4559	GHP Acrodisc GF, GF/0.45 µm, 25 mm	50/200
4472	Acrodisc CR PTFE, 0.45 µm, 4 mm	250/750
4552	Acrodisc CR PTFE Minispike, 0.2 µm, 13 mm	100/300
4423	Acrodisc CR PTFE, 0.2 µm, 13 mm	100/300
4553	Acrodisc CR PTFE Minispike, 0.45 µm, 13 mm	100/300
4422	Acrodisc CR PTFE, 0.45 µm, 13 mm	100/300
4225	Acrodisc CR PTFE, 0.2 µm, 25 mm	50/200
4219	Acrodisc CR PTFE, 0.45 µm, 25 mm	50/200
4226	Acrodisc CR PTFE, 1 µm, 25 mm	50/200
4484	Nylon Acrodisc, 0.45 µm, 4 mm	250/750
4550	Nylon Acrodisc Minispike, 0.2 µm, 13 mm	100/300
4427	Nylon Acrodisc, 0.2 µm, 13 mm	100/300
4551	Nylon Acrodisc Minispike, 0.45 µm, 13 mm	100/300
4426	Nylon Acrodisc, 0.45 µm, 13 mm	100/300
4436	Nylon Acrodisc, 0.2 µm, 25 mm	50/200
4438	Nylon Acrodisc, 0.45 µm, 25 mm	50/200
4549	Nylon Acrodisc GF, GF/0.45 µm, 25 mm	50/200
4450	Acrodisc LC PVDF Minispike, 0.2 µm, 13 mm	100/300
4455	Acrodisc LC PVDF, 0.2 µm, 13 mm	100/300
4452	Acrodisc LC PVDF Minispike, 0.45 µm, 13 mm	100/300
4457	Acrodisc LC PVDF, 0.45 µm, 13 mm	100/300
4406	Acrodisc LC PVDF, 0.2 µm, 25 mm	50/200
4408	Acrodisc LC PVDF, 0.45 µm, 25 mm	50/200

### Ion Chromatography Certification

Pall Life Sciences certifies the following Ion Chromatography (IC) Acrodisc syringe filters have been tested using a highly sensitive ion chromatography technique. To monitor inorganic extractables, statistically representative samples are tested. Actual background levels of filter extractables are typically less than 20 ppb for chloride, 6 ppb for nitrate, 1 ppb for phosphate, and 10 ppb for sulfate.

### Ordering Information

Part Number	Description	Pkg Pk/Case
4483	IC Acrodisc (PES), 0.2 µm, 13 mm	100/300
4485	IC Acrodisc (PES), 0.45 µm, 13 mm	100/300
4583	IC Acrodisc (PES), 0.2 µm, 25 mm	50/200
4585	IC Acrodisc (PES), 0.45 µm, 25 mm	50/200

### Syringe Filters for Prefiltration

Pall Life Sciences Glass Fiber Acrodisc syringe filters provide fast flow rates and high particulate holding capacities, which make it an excellent filter for prefiltration of heavily contaminated samples.

### Ordering Information

Part Number	Description	Pkg Pk/Case
4523	Glass Fiber Acrodisc, 1 µm, 25 mm	50/200
4524	Glass Fiber Acrodisc, 1 µm, 37 mm	15/60

### Instructions for Use

Before filling the syringe with sample, draw approximately 1 mL of air into the syringe. This will allow the air to follow the sample out of the syringe. This "air purge" minimizes fluid retention within the filter device.

Fill the syringe with the solution to be filtered. CAUTION: Use of syringes smaller than 10 mL can generate excessive pressure on the filter, which may exceed maximum operating pressure.

Holding the filter device in one hand and the filled syringe in the other, secure (without excessive force) the filled syringe to the filter device with a twisting motion.

Apply gentle pressure to begin filtration. (A gentle pressure helps assure maximum throughput.) CAUTION: As the filter removes particulate, filtration will become more difficult (the syringe plunger will be harder to use) and pressure will rapidly increase on the filter. Change filters when resistance becomes excessive. Failure to change filter may result in housing rupture, which results in particulate contaminating the filtrate.

These filters are for **single use only**.

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### R = RESISTANT

No significant change was observed in flow rate or bubble point.

### LR = LIMITED RESISTANCE

Moderate changes in physical properties of the membrane were observed.

### NR = NOT RESISTANT

The membrane is basically unstable and is not recommended for use.

### - = INSUFFICIENT DATA

Information not available. Trial testing is recommended.

‡ = UV absorbance was

set at 254 nm.

## Specifications

### Materials of Construction

Filter Media:

GHP: GH Polypro membrane (GHP, hydrophilic polypropylene)

PTFE: Hydrophobic polytetrafluoroethylene membrane on a polypropylene support

Nylon: Hydrophilic nylon membrane

PVDF: Hydrophilic polyvinylidene fluoride membrane

PES: Hydrophilic polyethersulfone Supor<sup>®</sup> membrane

Glass: Binder-free borosilicate glass fiber

Housings: Polypropylene

### Pore Size

0.2, 0.45, and 1 µm

### Effective Filtration Area

4 mm: 0.28 cm<sup>2</sup>

13 mm: 1.0 cm<sup>2</sup>

25 mm: 2.8 cm<sup>2</sup>

37 mm: 7.5 cm<sup>2</sup>

### Sample Volume

4 mm: < 2 mL

13 mm: < 10 mL

25 mm: < 100 mL

37 mm: < 200 mL

### Inlet/Outlet Connections

Female Luer-Lok<sup>®</sup> inlet/male luer outlet or minispike outlet as noted

### Typical Hold-up Volume (with air purge)

4 mm: < 10 µL

13 mm minispike: < 14 µL

13 mm: < 30 µL

25 mm: < 100 µL

25 mm GHP GF and Nylon GF: < 125 µL

37 mm: < 500 µL

### Maximum Operating Temperature

GHP, GHP GF, Nylon, Nylon GF, PES, and 13 mm PVDF: 55 °C (131 °F)

PTFE, and 25 mm PVDF: 100 °C (212 °F)

Glass Fiber: 135 °C (275 °F)

### Maximum Operating Pressure

PVDF: 3.4 bar (340 kPa, 50 psi)

Glass Fiber, PES, 4 mm Nylon, and 4 mm PTFE: 5.2 bar (520 kPa, 75 psi)

GHP GF, Nylon GF, and 25 mm Nylon: 6.2 bar (620 kPa, 90 psi)

GHP, PTFE, and 13 mm Nylon: 6.9 bar (690 kPa, 100 psi)

## WARNING

Employment of the products in applications not specified, or failure to follow all instructions contained in this product information insert, may result in improper functioning of the product, personal injury, or damage to property or the product. See Statement of Warranty in our most recent catalog.

The HPLC and Ion Chromatography certification does not apply to packaged products that have been opened or damaged by conditions outside of the control of Pall Life Sciences.

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Pall, , Acrodisc, and Supor are trademarks of Pall Corporation. © indicates a registered trademark in the USA. •Luer-Lok is a trademark of Becton Dickinson and Company.