### MicroFunnel™ Plus and MicroFunnel Plus AP Filter Funnels

- Disposable, ready-to-use sample container and filter funnel all in one.
- Vented lid allows unit to filter without removing lid, protecting sample from contamination.
- AP sampling port for collecting sample without opening lid.

#### Ordering Information

Prod. No.	Description	Packaging
<b>MicroFunnel Plus filter funnels</b>		
4807 <sup>†</sup>	MicroFunnel Plus filter funnel with 0.45 µm GN-6 Metrical® membrane, 100 mL capacity, individually bagged	50/pkg
4808 <sup>†</sup>	MicroFunnel Plus filter funnel with 0.45 µm Metrical Black membrane, 100 mL capacity, individually bagged	50/pkg
4809	MicroFunnel Plus filter funnel with 0.2 µm Supor® membrane, 100 mL capacity, individually bagged	50/pkg
4823	MicroFunnel Plus filter funnel with 0.45 µm Supor membrane, 100 mL capacity, individually bagged	50/pkg
4813	MicroFunnel 300 Plus filter funnel with 0.2 µm Supor membrane, 300 mL capacity, individually bagged	20/pkg
4814	MicroFunnel 300 Plus filter funnel with 0.45 µm Supor membrane, 300 mL capacity, individually bagged	20/pkg
4829 <sup>†</sup>	MicroFunnel 300 Plus filter funnel with 0.45 µm GN-6 Metrical membrane, 300 mL capacity, individually bagged	20/pkg

#### MicroFunnel Plus AP filter funnels

4843	MicroFunnel Plus AP with 0.2 µm Supor membrane, 100 mL capacity, individually bagged	50/pkg
4844	MicroFunnel Plus AP with 0.45 µm Supor membrane, 100 mL capacity, individually bagged	50/pkg

#### Accessories

4824	Funnel Holder, 100 mL	1/pkg
4825	Funnel Holder, 300 mL	1/pkg
4845	AP Sampling Tube, individually bagged, sterile	50/pkg

<sup>†</sup>PN 4807, 4808, and 4829 are not compatible with hot water loop applications.

#### Membranes

**GN-6 Metrical membrane:** A 0.45 µm mixed cellulose ester membrane used in the Membrane Filter (MF) Technique for the determination of Total Coliforms, *E. coli*, yeast, mold, Streptococcal strains, and total organism contamination.

**Metrical Black membrane:** A 0.45 µm modified polyethersulfone membrane used with the MF Technique to provide better contrast to light colored colonies of bacteria, yeast and mold.

**Supor 200 membrane:** A 0.2 µm polyethersulfone membrane used in the MF Technique for the determination of contamination by *Pseudomonas* species.

**Supor 450 membrane:** A 0.45 µm polyethersulfone membrane used in the MF Technique for the determination of total organism contamination.

The membrane contained in MicroFunnel filter funnels has been manufactured for microbiological analysis in accordance with the Membrane Filtration procedures referenced in "Standard Methods for the Examination of Water and Wastewater," 20th Edition, and the US Environmental Protection Agency's "Microbiological Methods for Monitoring the Environment," 600/8-78-017.

- Recovery testing has shown that the grid lines do not enhance or inhibit growth of organisms.
- The sterilization process does not enhance or inhibit subsequent growth of organisms.
- Recovery of total and fecal coliform for 0.45 µm GN-6 Metrical membrane, *Saccharomyces cerevisiae* for 0.45 µm Metrical Black membrane, *Pseudomonas aeruginosa* for 0.2 µm Supor membrane, and *Escherichia coli* for 0.45 µm Supor membrane indicates that there is no influence on microbial growth and development due to chemical extractables.

	GN-6 Metrical Membrane	Metrical Black Membrane	Supor 200 Membrane	Supor 450 Membrane
Membrane Diameter	47 mm	47 mm	47 mm	47 mm
Pore Size	0.45 µm	0.45 µm	0.2 µm	0.45 µm
Recovery (%) in comparison to a pour plate	≥ 90	≥ 85	≥ 90	≥ 90
Thickness	6 mils (152 µm)	5.1 mils (130 µm)	5.7 mils (145 µm)	5.5 mils (140 µm)
Typical Water Flow Rates	> 65 mL/min/cm <sup>2</sup> at 0.7 bar (70 kPa, 10 psi)	> 34	> 19	> 38

#### Specifications

##### Materials of Construction

Filter media: PN 4807, 4829 GN-6 Metrical (mixed cellulose ester) membrane, white  
PN 4808 Metrical Black (modified polyethersulfone) membrane, black  
PN 4809, 4823, 4813, 4814, 4843, 4844 Supor® (polyethersulfone) membrane, white

Support pad: Cellulose

Funnel, Base and Lid Plug: Polypropylene

Vented Lid: Polyethylene

Lid and Base Vent Membrane: hydrophobic Versapor® (acrylic copolymer on a non-woven support) membrane

Bag: Polyethylene

Adapter: Polyethylene

Adhesive Gasket: Urethane

##### Effective Filtration Area

13.46 cm<sup>2</sup>

##### Maximum Vacuum

635 mm Hg (25 in. Hg) (Vacuum Use Only)

##### Maximum Operating Temperature

4809, 4823, 4813, 4814: 90 °C (194 °F)

4807, 4808, 4829: Ambient

##### Sterilization

Sterilized by gamma irradiation

##### Dimensions

100 mL capacity:

Height: 7.6 cm (3.0 in.)

Diameter: 6.1 cm (2.4 in.)


300 mL capacity:

Height: 8.9 cm (3.5 in.)


Diameter: 8.7 cm (3.4 in.)

Figure 1

Prod. No. 4807, 4808, 4809, 4823, 4813, 4814 and 4829

Prod. No. 4843 and 4844

#### Instructions for Use

##### Collecting Ambient Samples, follow this sample collection procedure:

1. Select MicroFunnel™ Plus filter funnel based on temperature and volume of sample to be collected and pore size of membrane.
2. Remove funnel from the protective sealed bag by firmly grasping the bag above and below the tear notch. Evenly tear the bag by pulling up and over with a rotation of the wrist of the hand grasping above the tear notch. Avoid pulling the bag toward or away from you or trying to open the bag by pulling the sealed seam apart as this will cause the polyethylene plastic of the bag to stretch and tear unevenly.
3. a) MicroFunnel Plus filter funnel: Position the funnel near the sample point and remove the lid by placing the thumb under the tab on the lid. Be sure not to contaminate the funnel by touching the inner surfaces of the lid and funnel.
   
b) MicroFunnel Plus AP filter funnel: Aseptically open the adapter end of the sampling tube (PN 4845) and remove the lid plug from the sample port. Press the male luer adapter into the sample port. Remove the bag from the sampling tube and connect the end of the tube to the sanitized port on the water system.
4. a) MicroFunnel Plus filter funnel: Collect the sample and replace lid. Be sure the lid is firmly snapped down to the second liquid-tight engagement by applying downward pressure around the perimeter of the lid. This can be achieved by placing the funnel on a flat, sturdy surface and evenly pressing the lid down around the perimeter or grasping the funnel in a vertical position with thumbs on the upper perimeter of the lid and finger tips on the lower edge of the base and firmly squeezing the funnel until you hear the lid engage.
   
b) MicroFunnel Plus AP filter funnel: After the sample is collected, remove the sampling tube and replace the lid plug.
5. Return the filled funnels to the laboratory for filtration and analysis.

##### Collecting a Hot Water Sample, follow this sample collection procedure:

**NOTE: When using this product for sampling hot water loops, observe safety precautions which include the use of insulated rubber gloves, safety glasses, and Pall Life Sciences Funnel Holder for critical sample sites (PN 4824 for 100 mL funnels and PN 4825 for 300 mL funnels).**

1. Select MicroFunnel Plus filter funnel based on temperature, volume of sample to be collected, and pore size of membrane. Remove funnel from the protective sealed bag as described in step 2 above.
2. Place funnel into the holder with the graduations visible through the viewing slot.
3. a) MicroFunnel Plus filter funnel: Position the holder with the funnel near the sample point and remove the lid as described in step 3a above.
   
b) MicroFunnel Plus AP filter funnel: Aseptically open the adapter end of the sampling tube (PN 4845) and remove the lid plug from the sample port. Press the male luer adapter into the sample port. Remove the bag from the sampling tube and connect the end of the tube to the sanitized port on the water system.
4. a) MicroFunnel Plus filter funnel: Collect the sample and replace the lid. Be sure the lid is firmly snapped down to the second liquid-tight engagement as described in step 4a above. Sealing the lid is best accomplished by placing the holder, with the funnel inside the holder, on a flat sturdy surface and applying even pressure around the perimeter of the lid.
   
b) MicroFunnel Plus AP filter funnel: After the sample is collected, remove the sampling tube and replace the lid plug
5. Remove the funnel from the holder by pressing upward with your index finger through the opening in the bottom of the holder and gently lifting the funnel out. Be careful not to damage the vent tab on the bottom of the funnel. Return the filled funnels to the laboratory for filtration and analysis.

#### Filtration and Analysis

6. Prepare workspace and vacuum manifold set-up.
7. Position the funnel over the port of the manifold and remove the base vent membrane tab on the bottom of the funnel and place funnel directly on the manifold.
8. Turn on Vacuum to draw the sample through the membrane then turn off vacuum.
9. Remove lid and discard.
10. Remove funnel cylinder by grasping the base of the funnel and firmly squeezing the funnel cylinder at or below the midpoint of the cylinder. Squeeze the cylinder until you hear a slight snap of the adhesive gasket. Then tilt the cylinder either toward or away from you, lift the cylinder and discard.
11. Using aseptic technique, remove the membrane from the base.
12. Place membrane filter on a prepared agar plate or broth soaked pad in a Petri dish.
13. Invert and incubate.

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Pall, **PALL**, Metricel, MicroFunnel, Supor, Vacushield, and Versapor are trademarks of Pall Corporation. © Indicates a registered trademark in the USA.
© 2006, Pall Corporation, 06/06