AcroPrep™ 24 Filtration System

Description

The AcroPrep 24 filtration system is manufactured for use with the Waters Alliance® sample carousel. The AcroPrep 24 filtration system, consisting of a manifold and an AcroPrep disposable filter-plate, allows for the simultaneous filtration of 24 samples directly into an Alliance HPLC system sample carousel in only seconds.

Ordering Information

AcroPrep Filter-Plates

Part Number	Description	Pkg
186000154	GHP filter-plate, 0.45 µm	10/pkg
186000154P	Glass prefilter/GHP filter-plate, 0.45 µm	10/pkg
186000155	Nylon filter-plate, 0.45 µm	10/pkg
186000156	PTFE filter-plate, 0.45 µm	10/pkg
186000157	PVDF filter-plate, 0.45 µm	10/pkg
600000158	GHP filter-plate, 0.45 µm	2/pkg
186000158	GHP filter-plate, 0.2 µm	10/pkg
186000159	Nylon filter-plate, 0.2 µm	10/pkg
186000160	PTFE filter-plate, 0.2 µm	10/pkg

Manifold and Spare Parts

Part Number	Description	Pkg
289000159	Filtration manifold	1/pkg
700000231	Glass housing for filtration manifold	1/pkg
700000232	O-ring kit	3/pkg
700000233	Rubber feet	3/pkg
700000234	Manifold alignment post	3/pkg

Specifications

Materials of Construction

Housing: Polypropylene
Filter Media: GH Polypro (hydrophilic polypropylene), Nylon, PTFE, PVDF, or Glass Prefilter (borosilicate glass).

Pore Size

0.2 and 0.45 µm

Filter-plate Diameter

15.24 cm (6.0 in.) plate containing 24 wells

Effective Membrane Filtration Area per Well

0.94 cm²

Sample Well Volume

1.9 mL (recommended maximum sample volume: 1.7 mL)

Outlet Connections

Minispike outlet

Typical Hold-up Volume per Well

< 50 µL

Operating Temperature Range

18-27 °C (64-80 °F); Ambient

Maximum Vacuum

56.0 cm (22 in.) Hg

Figure 1



Figure 2

Instructions for Use

The AcroPrep 24 filtration system operates by vacuum. Conventional lab vacuum pumps such as Pall Life Sciences 115V vacuum pump (PN 13157) or 230V vacuum pump (PN 13158) are ideal for use. To obtain maximum convenience, we recommend the use of snap-cap vials to eliminate the need to remove vials for screw capping after filtration.

- Attach vacuum hose to the barbed connector on the base.
- Place your Alliance sample carousel with empty 2 mL vials into the manifold. Turn the sample carousel until it seats onto the center hub. The sample carousel is now properly positioned.
- Place the filter-plate onto the manifold matching the notches on the filter-plate with the alignment posts to orient the filter-plate over the sample vials in the carousel. The filter-plate has a numerically sequenced label which matches with the carousel's sample numbers. Make sure that the flange of the plate contacts the whole circumference of the O-ring.
- Apply vacuum to ensure that all the seals are working properly. Shut off vacuum source. Once vacuum is fully relieved, pipette samples into the filter wells of the AcroPrep filter-plate. Do not pipette more than 1.7 mL (1700 µL), otherwise sample could spill over into the adjacent filter well.
- Apply vacuum. Do not exceed a vacuum level of 56.0 cm (22 in.) Hg.
- Allow vacuum to run until all filter wells are evacuated. Continue running vacuum for 10-20 seconds, to eliminate potential dripping of samples from tips when the plate is removed.
- Shut off vacuum source. Once vacuum is relieved, lift the used filter-plate off the manifold and dispose of it.
- Cap sample vials and place the sample carousel into the Waters Alliance HPLC system or simply transfer the filtered sample to any HPLC sample carousel that houses 12 x 32 mm vials.

Instructions for replacing glass housing and O-rings:

To replace glass housing

Do not apply vacuum if glass is cracked or chipped, replace immediately.

- Unscrew and remove the three manifold alignment posts by turning them counterclockwise.
- Remove manifold support ring.
- Carefully remove glass housing.
- Place new glass housing (PN 700000231) onto manifold base plate so that the glass housing fits over the base O-ring.
- Place manifold support ring onto glass housing
- Insert alignment posts through holes in the support ring and hand-tighten them.

If your filtration system is not holding vacuum, you may have a damaged O-ring. Visually check O-rings for damage and replace if needed.

To replace top support O-ring

- Remove worn O-ring.
- Place the new O-ring (PN 700000232) into the O-ring groove.

Instructions for Use (continued)

To replace manifold base and support O-rings

- Unscrew and remove the three manifold alignment posts by turning them counterclockwise.
- Remove the manifold support ring.
- Carefully remove glass housing.
- Remove worn O-rings.
- Place the new O-rings (PN 700000232) into the O-ring grooves.
- Place glass housing onto manifold base plate so that glass housing fits over the base O-ring.
- Place manifold support ring onto glass housing using the alignment posts as a guide.
- Insert alignment posts through holes in the support ring and hand-tighten the three manifold alignment posts into the manifold base plate.

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Chemical Compatibility

Different polymeric membranes have different chemical compatibilities. Based on the application and chemical compatibility, there may be one or several AcroPrep filter-plate possibilities.

Solvent	GHP	Nylon	PVDF	PTFE
Acetone	R	R	NR	R
Acetonitrile	R	R	R	R
Acetic Acid, glacial	R	NR	R	R
n-Butanol	R	R	R	R
Chloroform	R	R	R	R
Dioxane	R	–	R	R
Dimethyl formamide	R	R	NR	R
Dimethyl sulfoxide	R	R	NR	R
Ethanol	R	R	R	R
Ethyl acetate	R	R	R	R
Ethyl ether	R	R	R	R
Freon TF	R	R	R	R
Hydrochloric acid (1N)	R	NR	R	R
Hexane, dry	R	R	R	R
Methanol	R	R	R	R
Methylene chloride	R	R	R	R
Methyl ethyl ketone	R	R	NR	R
N-Methylpyrrolidone	R	R	NR	R
Isopropanol	R	R	R	R
Sodium hydroxide (3N)	R	LR	NR	R
Tetrahydrofuran	R	R	NR	R
Tetrahydrofuran/Water (50/50)	R	R	–	R
Toluene	R	R	R	R
Water	R	R	R	R

R = RESISTANT **NR = NOT RESISTANT**
LR = LIMITED RESISTANT **-- = INSUFFICIENT DATA**

Test Methods: The data presented in this chart is a compilation of testing by Pall Life Sciences with certain chemicals, manufacturer's data, or compatibility recommendations from the Compass Corrosion Guide, by Kenneth M. Pruett. This data is intended to provide expected results when filtration devices are exposed to chemicals under static conditions for less than 48 hours at 25 °C (77 °F) unless otherwise noted. This chart is intended only as a guide. Accuracy cannot be guaranteed. Users should verify chemical compatibility with a specific filter under actual use conditions. Chemical compatibility with a specific filter, under actual use conditions, is affected by many variables, including temperature, pressure, concentration, purity, and various chemical combinations which prevent complete accuracy.

Complementary Products

To meet all your requirements for analytical sample filtration, Pall Life Sciences also offers its Acrodisc® syringe filter line with PTFE, GH Polypro, Nylon, and PVDF membranes. These filters are available in 25 mm with standard male luer outlets and in 13 mm with either male luer or minispike outlets.

Pall Life Sciences also offers 47 and 50 mm membrane disc filters for the purification and degassing of the mobile phase solvents used in liquid chromatography applications. The membrane disc filters are identical in composition and quality to those used in HPLC-certified AcroPrep filter-plates and Acrodisc syringe filters.

For additional information, contact your local laboratory distributor or Pall Life Sciences Technical Service Department at 800-521-1520 (USA), or 734-665-0651 (other).

WARNING

Employment of the products in applications not specified, or failure to follow all instructions contained in this product information insert, may result in improper functioning of the product, personal injury, or damage to property or the product. See Statement of Warranty in our most recent catalog.

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