



AcroSep™

Chromatography Columns

Act Fast! 40% discount for limited time only!

Accelerating protein purification and analysis

Purchase select AcroSep Chromatography columns and receive 40% off!

Hurry! Offer valid January 1 thru March 31, 2010 – order today!

AcroSep Chromatography Columns

- ▶ **Flexible.** More separation and elution options in a single format to provide greater purity.
- ▶ **Versatile.** Luer lock inlet and outlet allows convenient use with syringe, pump, or automated chromatography system.
- ▶ **User-friendly design.** Color-coded and labeled by chemistry type. Collar is hexagonal so columns will not unexpectedly roll off lab surface.

Powerful Performance Advantages

- ▶ **Protein A Ceramic HyperD® resin** offers high selectivity with low non-specific binding.
- ▶ Achieve high purity and yield of target proteins with flexible, uncharged **IMAC HyperCel™ resin**.
- ▶ **Blue Trisacryl® M resin** allows good diffusion and improves exchange kinetics.
- ▶ **MEP HyperCel resin** offers broad species and isotype binding capabilities for antibodies.
- ▶ Receive effective removal of a wide range of ionic and non-ionic detergents with **SDR HyperD resin**.

Qualifying Products

VWR Cat. No.	Description	*Pkg. of 5
MP97011-234	Protein A Ceramic HyperD F, pearl	\$252.00
MP97011-242	IMAC HyperCel, cobalt blue	\$ 66.15
MP97015-750	Blue Trisacryl M, dark blue	\$ 78.00
MP97011-236	MEP HyperCel, purple	\$ 88.20
MP97006-690	SDR HyperD, natural	\$ 96.00

*Price reflects 40% discount.

To order contact:



vwr.com
1.800.932.5000

1/10, GN09.3296POD



Prices, product appearance and specifications are current at the time of printing, subject to change without notice. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in USA only, void where prohibited by law or company policy, while supplies last. Visit vwr.com to view our privacy policy and additional disclaimers.

VWR, forms of VWR and the VWR logo and/or design are either registered trademarks ®, trademarks ™, or service marks ℠ of VWR International, LLC in the United States and/or other countries. All other marks referenced herein are registered trademarks, trademarks or service marks of their respective owner(s). For a complete list of trademark owners please visit vwr.com. ©2010 VWR International, LLC. All rights reserved.