



The only syringe filter

certified

for LCMS applications

Acrodisc® MS Syringe Filters

Designed specifically for liquid chromatography mass spectrometry (LCMS)

- ▶ **LCMS certified** – Minimize interference in your results. First LCMS filter certified for extremely low levels of extractables.
- ▶ **Low ion suppression/enhancement** – Reduce the need for retesting. Acrodisc MS syringe filters will not contribute extractables that interfere with the ionization process – the heart of the LCMS technique.
- ▶ **Protective packaging design** – Packaged into five separate tubes to maintain purity until ready to use.

Request your free sample today.

Offer valid: September 1–February 28, 2012

Product Ordering Information

VWR Cat. No.	Pall Part Number	Description	Pkg	List Price
89194-420	MS-3201	0.2 µm, 25 mm, WWPTFE membrane, LCMS certified	50/pkg	\$183.50

To receive your FREE sample, visit www.vwrsp.com/promotions, search for Promo Code 3315 and fill in the online form by February 28, 2012.



Prices, product appearance, and specifications are current at the time of printing, subject to change without notice. Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in U.S. dollars unless otherwise noted. Offers valid in USA only, void where prohibited by law or company policy, while supplies last. Visit vwr.com to view our privacy policy and additional disclaimers.

VWR, forms of VWR, and the VWR logo and/or design are either registered trademarks ®, trademarks ™, or service marks SM of VWR International, LLC in the United States and/or other countries. All other marks referenced herein are registered trademarks, trademarks, or service marks of their respective owner(s). For a complete list of trademark owners please visit vwr.com.

©2011 VWR International, LLC. All rights reserved.