



Medical

400 mL Transfer Bag with Multileads – Designed for Cryoprecipitate Pooling



A perfect fit for the Leukotrap® Synergy Platform systems package to help your center Improve Supply, Increase Safety, and Grow Savings.

- ▶ Six (6) manifold leads allow Blood Establishments to determine their own pool sizes.
- ▶ Pall's Transfer Bag offers a time saving, readily-available alternative to cryoprecipitate pooling on demand.
- ▶ All tubing is compatible with standard sterile tubing connection devices.
- ▶ Set design allows for pre-storage pooling to be performed in a "closed system".
- ▶ The 400 mL Transfer Bag can be packaged in standard FFP cartons which can easily be stored in freezer units and distributed in shipping containers with minimal breakage.
- ▶ Provides an opportunity to expand product mix and services.



400 mL Transfer Bag with Multileads

Features & Benefits

- ▶ **Spiked leads** allow flexibility for open or closed system processing.
- ▶ **7 inch manifold** reduces product loss and maximizes Cryoprecipitated AHF recovery.
- ▶ **6-to-1 manifold tubing connector** eliminates product loss associated with standard Y-connected tubing leads.
- ▶ **Temporary clamp** located between the manifold and transfer bag helps to control the pooling process.
- ▶ **400 mL Transfer Bag** is designed for storage in standard FFP containers.
- ▶ **Interim ISBT128 & Codabar label on enhanced paper** for improved adhesion of overlables.

Ordering Information

- ▶ **Part Number:** 732-42M
- ▶ **Description:** 400 mL Transfer Bag with 6 Multileads
- ▶ **Packaging:** 96 sets, bulk packaged (24 sets/inner box; 4 inner boxes per case)

Latex content: This product is free of natural rubber latex.

Contact your **Pall Medical Account Manager** or call our **Customer Care Hotline at 800.645.6578** to obtain more information about the 400 mL Transfer Bag with Multileads.



United States
25 Harbor Park Drive
Port Washington, NY 11050

1.800.645.6578 toll free phone (USA)
1.516.484.5400 phone
1.516.801.9754 fax




Visit us on the Web at www.pall.com/cryopooling

International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States and Venezuela. Distributors in all major industrial areas of the world.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

© 2010, Pall Corporation. Pall, , and Leukotrap are trademarks of Pall Corporation. ® indicates a trademark registered in the USA. **Redefining the Business of BloodSM** is a service mark of Pall Corporation.