



Medical



Pall PharmAssure Pharmacy Filtration Devices

25 mm Hydrophobic Vent Filter

| | |
|--------------------------------------|---|
| Housing | Acrylic |
| Membrane | 0.2 µm Hydrophobic acrylic co-polymer with polyamide support media |
| Effective Filtration Area | 2.8 cm ² |
| Connections | ISO Female Luer Lock Inlet ISO Male Luer Lock Outlet |
| Maximum Operating Pressure | 15 psi (1 bar) |
| Maximum Operating Temperature | 55°C |
| Endotoxin Level | <0.25 EU/mL (Limulus Amoebocyte Lysate (LAL) test) |
| Biocompatibility | Components meet USP Biological Reactivity Testing <i>in vivo</i> <88> |
| Sterilisation Method | Gamma irradiation |

Reorder Information

| Product Code | Pore Size | Colour | Packaging |
|--------------|-----------|-------------|-----------|
| HP2002 | 0.2 µm | Clear/Clear | 50/pkg |

Technical Data Sheet



Medical

Europa House, Havant Street
Portsmouth PO1 3PD, England

+44 (0)23 9230 3452 telephone
+44 (0)23 9230 3541 fax
Biosvc@Pall.com E-mail



Visit us on the web at www.pall.com

International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, Hong Kong, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States and Venezuela. Distributors in all major industrial areas of the world.

This document is not for distribution in the USA and Canada.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly. Part numbers quoted above are protected by the Copyright of Pall Europe Ltd.

Pall and PharmAssure are trade marks of Pall Corporation
Filtration. Separation. Solution. is a service mark of Pall Corporation
©2000, 2001 Pall Europe Limited

Printed in England. PMED/2M/DBD/0203

Filtration. Separation. Solution.SM



Medical

Technical Data Sheet

Pall PharmAssure Pharmacy Filtration Devices

Dispensing Pin



| | |
|--------------------------------------|--|
| Housing | Acrylic |
| Vent Membrane | 0.2 µm PTFE (polytetrafluoroethylene) |
| Protective Caps | Polyethylene |
| Effective Vent Area | 1.0 cm ² |
| Connections | ISO Female Luer Lock Inlet Dual Lumen Male Spike |
| Water Breakthrough | 5 psi (0.35 bar) |
| Maximum Operating Temperature | 55°C |
| Endotoxin Level | <0.25 EU/mL (Limulus Amoebocyte Lysate (LAL) test) |
| Biocompatibility | Components in fluid pathway meet USP Biological Reactivity Testing <i>in vivo</i> <88> |
| Sterilisation Method | Gamma irradiation |

Reorder Information

| Product Code | Vent Size | Colour | Packaging |
|--------------|-----------|--------|-----------|
| HP2402 | 0.2 µm | Opaque | 50/pkg |



Medical

Europa House, Havant Street
Portsmouth PO1 3PD, England

+44 (0)23 9230 3452 telephone
+44 (0)23 9230 3541 fax
Biosvc@Pall.com E-mail



Visit us on the web at www.pall.com

International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, Hong Kong, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States and Venezuela. Distributors in all major industrial areas of the world.

This document is not for distribution in the USA and Canada.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly. Part numbers quoted above are protected by the Copyright of Pall Europe Ltd.

Filtration. Separation. Solution.SM

Pall and PharmAssure are trade marks of Pall Corporation
Filtration. Separation. Solution. is a service mark of Pall Corporation
©2000, 2001 Pall Europe Limited

Printed in England. PMED/2M/DBD/0203