

DECLARATION OF COMPLIANCE

Seitz® HR Series Filter Sheets

Product Grades: HR100, HR2000 and HR9000

The Seitz HR Series filter sheets are comprised of a cellulose matrix with perlite and diatomaceous earth. They also contain binder resin and polyolefin fibers.

These filter sheets are intended for filtration of aqueous solutions including cellulase enzymes. Seitz HR filter sheets may be used for non-alcoholic as well as alcoholic beverages.

An initial flush is recommended prior to use.

Issued	1 September 2009
Revised	4 October 2024
Expires	28 February 2027
Reference	FBDCHRENm
Page	1 of 2



Deniz Karsu
Quality Assurance & Regulatory Affairs Manager
Pall Filtersystems GmbH

Components

Cellulose, binder resin, perlite, diatomaceous earth and polyolefin fibers

Declaration

Europe

Seitz HR Series filter sheets meet the requirements of food contact as detailed in European Regulation (EC) Number 1935/2004 in that:

- Our suppliers state that the monomers and additives of the polymer fibers are listed in European Regulation (EU) Number 10/2011 Annex I.
- The product complies with German Recommendations XXXVI and XXXVI/1, as well as with the German Foodstuffs and Animal Feed Code (LFGB) §§ 30 and 31.

Users should satisfy themselves that these materials are suitable for use in their specific food application.

Mercosur

The Seitz HR Series filter sheets meet the requirements of food contact as detailed in Reglamento Técnico Mercosur sobre materiales celulósicos para cocción y filtración en caliente, Mercosur/GMC/Res. No. 41/15.

USA

The following materials of construction meet the FDA requirements for food contact use as detailed in the US Code of Federal Regulations, 21 CFR paragraphs 177-199 in that:

- Cellulose and binder resin to 21 CFR section 177.2260 (Filters, resin bonded) and 21 CFR section 176.170 (Components of paper and paperboard in contact with aqueous and fatty foods)
- Polyolefin fibers to 21 CFR section 177.1520 (Olefin polymers)
- Total extractables as per 21 CFR 177.2260 (Filters, Resin Bonded) (g) (h) (i) (j) (k) (l).
50% ethanol at room temperature and n-hexane at reflux were used in the extractables testing.

The following are listed in the Food Chemical Codex (FCC): Perlite and diatomaceous earth

Process Quality System

Site of Manufacture: Pall Filtersystems GmbH, Bad Kreuznach, Germany supplied by Pall Sàrl. Made in Germany.

The Quality Management System at Pall Filtersystems GmbH, Bad Kreuznach is certified to ISO 9001:2015.

These products / product packaging carry a lot number / date code to facilitate traceability to suppliers' materials and Pall production records.

Pall Filtersystems GmbH confirm that this product is manufactured in line with the principles of food contact materials GMP as detailed in Regulation 2023/2006.

Supplied in Europe by

Pall Sàrl
Av. de Tivoli 3
Fribourg
Switzerland
CH-1700

Visit us on the Web at www.pall.com/foodandbev



+1-866-905-7255 **Food and Beverage toll free**
foodandbeverage@pall.com

Corporate Headquarters

Port Washington, NY, USA
+1-800-717-7255 toll free (USA)
+1-516-484-5400 phone

European Headquarters

Fribourg, Switzerland
+41 (0)26 350 53 00 phone


Asia-Pacific Headquarters

Singapore
+65 6389 6500 phone

Pall Corporation has offices and plants throughout the world. To locate the Pall office or distributor nearest you, visit www.pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

IF APPLICABLE Please contact Pall Corporation to verify that the product conforms to your national legislation and/or regional regulatory requirements for water and food contact use.

© Copyright 2024, Pall Corporation. Pall,  and Seitz are trademarks of Pall Corporation. ® Indicates a trademark registered in the USA.

FBDCHRENm
OCTOBER 2024