GeneDisc® and Microbiology Products

Accelerating microbiological detection with accurate, flexible methods







DAY	11	: 2	:3	:4	: 5	¦6	:7	¦ 8	:9	:10	
	1	1 I I I I I I I I I I I I I I I I I I I	1	1	1	1	1	1	1	1	I Product
	1	1 I I I I I I I I I I I I I I I I I I I	1	1	1	1	1	1	1	1	l Volume



APPLICATION WORKFLOW



Reduced Risk of Cross-Contamination and Improved Flexibility





Fits in a standard benchtop autoclave

Modular 3 or 6 place operation

··· Front loading valves





Configurable with various filter funnels





Improved Microbial Recovery on Pall Membrane vs. Alternative Membrane









Pall Metricel® Black PES Membrane Day 7 Average Colony Forming Unit 123 CFU

Alternative Supplier Nitrocellulose Membrane Day 7 Average Colony Forming Unit 91 CFU