

STATEMENT OF COMPLIANCE: ISO 12500-3:2009 Emflon® PFRW Filter Cartridges SBF Junior Style "W" Code

Cartridge Part Number

SBF1PFRW _____

This is a guide to the part numbering structure only. For specific options, please contact Pall.

Table 1 : Seal Material

| Code | Description |
|------|-------------|
| H4 | Silicone |
| | elastomer |
| J | Ethylene |
| | propylene |
| | rubber |

Emflon PFRW SBF Junior Style final air filters are designed for gas filtration in the food and beverage industry, fermenter inlet air and exhaust venting, sterile process air and sterile venting of tanks in food and beverage processing.

'Sterile air' is defined as that filtered using a $0.2~\mu m$ liquid rated filter (capable of removing 10^7 of Brevundimonas diminuta per cm² of effective filtration area, in laboratory tests). Please note this product is not supplied pre-sterilized.

The key feature of this product is the inclusion of two layers of hydrophobic, 0.2 µm microbial rated Emflon PTFE medium in a polypropylene hardware filter construction.

Issued 6 March 2013 Revised 11 April 2017 Expires 28 February 2019

Page 1 of 2

Steven Bailey Quality Manager

Pall Manufacturing (UK) Ltd., Ilfracombe

Emflon PFRW Filter Cartridges (SBF Junior Style "W" Code)

Components

Filter Media Two layers of hydrophobic Emflon PTFE media

Support / Drainage Polypropylene

Hardware Polypropylene

O-ring Seal Silicone Elastomer (H4)

Ethylene Propylene Rubber (J)

Declaration

Emflon PFRW filters meet the requirements of ISO12500-3:2009 in that:

Emflon PFRW filter cartridges (tested as part number: AB1PFR7WH4) have demonstrated applicable performance capability in laboratory tests as described in report LrF 03b/12 appended. These studies were conducted on randomly sampled product made with the same Pall proprietary hydrophobic PTFE filter membrane, polypropylene support /drainage material and polypropylene hardware; and using the same manufacturing techniques.

Supplied in Europe by

Pall International Sàrl Av. de Tivoli 3 Fribourg Switzerland CH-1700



Pall Food and Beverage

New York - USA +1 516 484 3600 telephone +1 866 905 7255 toll free foodandbeverage@pall.com

Visit us on the web at www.pall.com

Pall Corporation has offices throughout the world. For Pall representatives in your area, please go to www.pall.com/contact.

Because of developments related to products, systems and/or services described herein, the data and procedures are subject to change without notice. Please consult your Pall representative or visit www.pall.com to verify that this information remains valid.

© Copyright 2017, Pall Corporation. Pall, and Emflon are trademarks of Pall Corporation. ® indicates a trademark registered in the USA. Filtration.Separation.Solution is a service mark of Pall Corporation.