

STATEMENT OF COMPLIANCE: ISO 12500-3:2009
Emflon® PFRW Filter Cartridges
SBF Junior Style “W” Code

Cartridge Part Number

SBF1PFRW
Table 1

This is a guide to the part numbering structure only. For specific options, please contact Pall.

Table 1 : Seal Material

Code	Description
H4	Silicone elastomer
J	Ethylene propylene rubber

Emflon PFRW SBF Junior Style final air filters are designed for gas filtration in the food and beverage industry, fermenter inlet air and exhaust venting, sterile process air and sterile venting of tanks in food and beverage processing.

The key feature of this product is the inclusion of two layers of hydrophobic, 0.2 µm microbial rated Emflon PTFE medium in a polypropylene hardware filter construction.

Issued 6 March 2013
Revised 11 February 2021
Expires 28 February 2023
Page 1 of 2



Andrew Betts
Quality Manager
Pall Manufacturing (UK) Ltd., Ilfracombe

Emflon PFRW Filter Cartridges (SBF Junior Style “W” Code)

Components

Filter Media	Two layers of hydrophobic Emflon PTFE media
Support / Drainage	Polypropylene
Hardware	Polypropylene
O-ring Seal	Silicone Elastomer (H4) Ethylene Propylene Rubber (J)

Declaration

Emflon PFRW filters meet the requirements of ISO12500-3:2009 in that:

Emflon PFRW filter cartridges (tested as part number: AB1PFR7WH4) have demonstrated applicable performance capability in laboratory tests as described in report UN2-170524-TFF96800-041c. These studies were conducted on randomly sampled product made with the same Pall proprietary hydrophobic PTFE filter membrane, polypropylene support /drainage material and polypropylene hardware; and using the same manufacturing techniques.

Supplied in Europe by

Pall International Sàrl
Av. de Tivoli 3
Fribourg
Switzerland
CH-1700




Pall Food and Beverage

New York - USA
+1 516 484 3600 telephone
+1 866 905 7255 toll free
foodandbeverage@pall.com

Visit us on the web at www.pall.com

Pall Corporation has offices throughout the world. For Pall representatives in your area, please go to www.pall.com/contact.

Because of developments related to products, systems and/or services described herein, the data and procedures are subject to change without notice. Please consult your Pall representative or visit www.pall.com to verify that this information remains valid.

© Copyright 2021, Pall Corporation. Pall,  and Emflon are trademarks of Pall Corporation. ® indicates a trademark registered in the USA.