

DECLARATION OF COMPLIANCE

Profile® II Filter Cartridges

RMF Style “W” Code

Cartridge Part Number

RM F W

Table 1 Table 2

This is a guide to the part numbering structure only. For specific options, please contact Pall.

Table 1 : Nominal Length

Code	Description
1	25.4 mm (10")
2	508 mm (20")
29	736.6 mm (29")
3	762 mm (30")
4	1016 mm (40")

Table 2 : Removal Rating

Code	Description
005	0.5 µm
010	1 µm
030	3 µm
050	5 µm
100	10 µm
200	20 µm
400	40 µm

RMF Profile II filter cartridges comprise of a polypropylene filter medium cylindrical sleeve formed by a proprietary Pall process with a polypropylene molded core with integral end cap / seal closures.

Profile II RMF style filter medium does not contain surfactants or binder resins.

RMF Profile II filter cartridges are suitable for use in aqueous-based food and beverage applications, excluding fatty food.

An initial flush is recommended prior to use.

Issued 1 May 2009
 Revised 24 April 2025
 Expires 28 February 2027
 Reference FBDCPROFRMFENh
 Page 1 of 3



Syed Sirajuddin
Pall Sàrl and MEA QARA Manager

Profile II Filter Cartridges (RMF Style “W” Code)

Components

Filter Media	Polypropylene
Core	Polypropylene
Elastomeric End Cap / Seal	Thermoplastic Elastomer material

Declaration

Profile II RMF Style “W” Code filter cartridges comprise of materials that meet regulatory and legislative requirements and guidelines for food contact in that:

Europe

The “W” Code RMF Profile II filter cartridges meet the requirements for food contact as detailed in European Regulation (EC) Number 1935/2004 in that:

Our suppliers' state that the polymeric components of construction are made from monomers and additives consistent with Annex I and II of Commission Regulation (EU) Number 10/2011 and its amendments, relating to plastic materials and articles intended to come into contact with foodstuffs (excluding seals).

OML and SML migration testing was performed on Profile II filter medium and polypropylene hardware, as AB1A0507WH4 (excluding elastomeric endcap material), and met migration criteria after flushing and in flow conditions, in;

- Simulant B (3% acetic acid) for 2 hours at reflux
- Simulant D1 (65% ethanol) for 2 hours at reflux
- Distilled water for 2 hours at 70 °C (158 °F)

OML and SML migration testing of the elastomeric endcap material has been performed, and met migration criteria after flushing and in flow conditions, in;

- Simulant A (10% ethanol), up to 70 °C (158 °F), for 2 hours
- Simulant B (3% acetic acid), up to 70 °C (158 °F), for 2 hours
- Simulant C (20% ethanol), up to 70 °C (158 °F), for 2 hours

Note:

This product contains materials that are subject to Specific Migration Limit (SML) requirements.

This product contains calcium stearate and amorphous silicon dioxide, which is approved as a direct food additive.

Users should satisfy themselves that these materials are suitable for use in their specific food application.

Mercosur

The materials of construction meet the requirements for food contact as detailed in Reglamento Técnico Mercosur

sobre materiales lista positiva de

- monómeros, otras sustancias de partida y polímeros autorizado para la elaboración de envases y equipamientos plásticos in contacto con alimentos, Mercosur/GMC/Res. No. 19/21 y
- aditivos para materiales plásticos destinados a la elaboración de envases y equipamientos in contacto con alimentos, Mercosur/GMC/Res. No 56/92.

USA

The materials of construction meet the FDA requirements for food contact use as detailed in Code of Federal Regulations, 21 CFR paragraphs 170-199, in that:

- Polypropylene to 21 CFR section 177.1520 (Olefin polymers)
- Thermoplastic Rubber - Supplier states that the material is acceptable for use in food contact applications under 21 CFR §177.1810(B)(3) and Food Contact Notification (FCN) 679

Process Quality System

Manufactured in the UK and supplied by Pall Sàrl. Made in UK.

The Quality Management System at the manufacturing site is certified to ISO 9001:2015.

These products / product packaging carry a lot number / date code to facilitate traceability to suppliers' materials and Pall production records.

Pall Sàrl confirms that the product is manufactured in line with the principles of food contact materials GMP as detailed in Regulation 2023/2006.

Supplied in Europe by

Pall Sàrl
Av. de Tivoli 3
Fribourg
Switzerland
CH-1700



+1-866-905-7255 **Food and Beverage toll free**
foodandbeverage@pall.com

Corporate Headquarters

Port Washington, NY, USA
+1-800-717-7255 toll free (USA)
+1-516-484-5400 phone

European Headquarters

Fribourg, Switzerland
+41 (0)26 350 53 00 phone

Asia-Pacific Headquarters


Singapore
+65 6389 6500 phone

Visit us on the Web at www.pall.com/foodandbev

Pall Corporation has offices and plants throughout the world. To locate the Pall office or distributor nearest you, visit www.pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

IF APPLICABLE Please contact Pall Corporation to verify that the product conforms to your national legislation and/or regional regulatory requirements for water and food contact use..

© Copyright 2025, Pall Corporation. Pall,  and Profile are trademarks of Pall Corporation. ® Indicates a trademark registered in the USA.

FBDCPROFRMFENh
APRIL 2025