

FOR A LIMITED TIME

PURCHASE \$2000 of products and **GET UP TO \$1000** additional products for **FREE!***

Pall Laboratory's Analytical Product family is the choice for fast flow rates and high throughput

HPLC / UHPLC SAMPLE FILTRATION

Nanosep® Centrifugal Devices

- Spin multiple samples at one time
- For small sample applications
- Available in 0.2 µm for UHPLC and 0.45 µm for HPLC applications




Acrodisc® Syringe Filters

- For processing one sample at a time
- Available in 13 mm and 25 mm formats
- Wide range of membrane types (wwPTFE, Nylon, PVDF, PTFE, Glass Fiber, etc)
- Available with prefilter for difficult to filter samples
- Available in 0.2 µm (UHPLC) and 0.45 µm (HPLC)



AcroPrep™ Advance Filter Plates


- For high throughput processing
- Available in 350 µL, 1 mL, and 2 mL formats
- Available in 0.2 µm for UHPLC and 0.45 µm for HPLC applications



MOBILE PHASE FILTRATION



ENSURE EXCELLENT DATA QUALITY!



* Purchase a minimum of \$2000 of Pall Microbiology QC products in one order and get up to \$1000 of FREE product. Limit one redemption per customer. \$2000 purchase and \$1000 credit must be for different products. Value of free products is based on Pall list price. Products Ordered or Free must be from list associated for this promotion. Promotion ends December 31, 2022.

For further details or to redeem for your free product, contact your Pall Sales representative.

FEATURED PRODUCTS

Critical decisions are often based on the results of your analytical QC testing. Confidence in those results comes from using products that are proven to perform as expected.

Acrodisc Syringe Filters for your HPLC and UHPLC sample prep workflows. 0.45 µm (HPLC) and 0.2 µm (UHPLC) can extend column life up to 52x. 0.2 µm pore size for UHPLC can extend column life up to 117x.

Solvac® Filter Holder and membranes for mobile phase filtration. Quickly degas mobile phase for smooth baseline while eliminating the safety hazards of using a glass funnel.

Nanosep Centrifugal Devices with wwPTFE membrane for small volumes. Prepare samples of < 1 mL for HPLC / UHPLC.

AcroPrep Advance Advance Filter Plates for high throughput testing. Compatible with automated liquid handlers for high volume laboratories.

Top-Selling Products for Analytical Chemistry

Purchase \$2000 Worth of These Products and get up to an Additional \$1000 Worth for Free!

Syringe Filters for UHPLC*

Part Number	Description	Pkg
AP-4910	Acrodisc One wwPTFE 25 mm 0.2 µm	50/pkg, 200/cs
AP-4913	Acrodisc One GXF wwPTFE 25 mm 0.2 µm	50/pkg, 200/cs
2400	Acrodisc One wwPTFE 13 mm 0.2 µm	100/pkg, 300/cs

High Throughput for UHPLC*

Part Number	Description	Pkg
8582	AcroPrep Advance wwPTFE 350 µL 0.2 µm	10/pkg
8682	AcroPrep Advance wwPTFE 1 mL 0.2 µm	5/pkg
8782	AcroPrep Advance wwPTFE 2 mL 0.2 µm	5/pkg

Small Volume for UHPLC

Part Number	Description	Pkg
ODPTFE02C34	Nanosep MF wwPTFE 0.2 µm	100/pkg

Syringe Filters for HPLC*

Part Number	Description	Pkg
AP-4916	Acrodisc One wwPTFE 25 mm 0.45 µm	50/pkg, 200/cs
AP-4919	Acrodisc One GXF wwPTFE 25 mm 0.45 µm	50/pkg, 200/cs
2402	Acrodisc One wwPTFE 13 mm 0.45 µm	100/pkg, 300/cs

High Throughput for HPLC*

Part Number	Description	Pkg
8584	AcroPrep Advance wwPTFE 350 µL 0.45 µm	10/pkg
8684	AcroPrep Advance wwPTFE 1 mL 0.45 µm	5/pkg
8784	AcroPrep Advance wwPTFE 2 mL 0.45 µm	5/pkg

Small Volume for HPLC

Part Number	Description	Pkg
ODPTFE04C34	Nanosep MF wwPTFE 0.45 µm	100/pkg

Syringe Filters for MS

Part Number	Description	Pkg
MS-3201	Acrodisc MS wwPTFE 25 mm 0.2 µm	50/pkg
MS-3301	Acrodisc MS wwPTFE 13 mm 0.2 µm	60/pkg

Syringe Filters for IC

Part Number	Description	Pkg
4583	Acrodisc IC Supor 25 mm 0.2 µm	50/pkg, 200/cs

Filters for Mobile Phase Filtration

Solvac Filter Holder

Part Number	Description	Pkg
4020	Solvac Filter Holder	1/pkg

UHPLC*

Part Number	Description	Pkg
66602	Nylon 47 mm 0.2 µm	100/pkg
60539	wwPTFE 47 mm 0.2 µm	50/pkg
66143	PTFE 47 mm 0.2 µm	100/pkg

HPLC*

Part Number	Description	Pkg
66609	Nylon 47 mm 0.45 µm	100/pkg
60548	wwPTFE 47 mm 0.45 µm	50/pkg
66149	PTFE 47 mm 0.45 µm	100/pkg

*other medias available

Includes any Acrodisc product with a list price below \$1000



PALL CORPORATION

Corporate Headquarters

25 Harbor Park Drive
Port Washington, New York 11050

Visit us on the Web at www.pall.com/lab

E-mail us at LabCustomerSupport@pall.com

© 2022 Pall Corporation. Pall, PALL, Acrodisc, AcroPrep, Nanosep, and Solvac are trademarks of Pall Corporation. ® indicates a trademark registered in the USA. Limit one per customer. Offer not valid where prohibited by law or institutional policy. Valid only in both US and CA. Allow up to 8 weeks for delivery of item. Offer expires December 31, 2022.