

Ordering Information

AB			Z		7				
Code	Nominal Length	Code	Removal Rating	Filter Area ⁽⁸⁾	Cartridge Style ⁽¹⁰⁾	Code	Filter Grade	Code	Gasket Option
1	10 in. (254 mm)	U010	1 µm	0.50 m ² (5.4 ft ²)	Double 226 O-ring with bayonet lock and fin end	P	Code outside US	H4	Silicone
2	20 in. (508 mm)	U2-20 ⁽⁹⁾	2 µm	0.65 m ² (7.0 ft ²)		Omit	Code in US only		
3	30 in. (762 mm)	U030	3 µm	0.56 m ² (6.0 ft ²)		Supplied with Certificate of Test and Full Traceability			
4	40 in. (1016 mm)	U6-40 ⁽⁹⁾	6 µm	0.56 m ² (6.0 ft ²)					
		U100	10 µm	0.49 m ² (5.3 ft ²)					
		U200	20 µm	0.40 m ² (4.4 ft ²)					
		U400	40 µm	0.39 m ² (4.2 ft ²)					

⁽⁸⁾ Nominal filter area per 10 in. (254 mm) filter element.

⁽⁹⁾ Multi-layered graded composites; 20 to 2 µm and 40 to 6 µm respectively.

Z in part number designates positive zeta in aqueous service.

⁽¹⁰⁾ Other adapter codes available.

Other materials available on request.